

**SPONTANEOUS CORRECTION OF
AN ALVEOLAR DEFORMATION**

Dr. Patrick Fellus



**THIS PATIENT ARRIVED AT THE AGE OF 5 (PHOTO TOP LEFT) WITH
DYSFUNCTION ASSOCIATED WITH ALVEOLAR DEFORMATION.**

**IN TWO MONTHS (PHOTO TOP RIGHT) WE RESTORED A FUNCTIONAL BALANCE
THANKS TO WEARING THE FROGGY MOUTH ALONE AND THE DEFORMATION
WAS REDUCED BY MORE THAN HALF.**



HE THEN EXPERIENCED A SPONTANEOUS EVOLUTION BETWEEN HIS 5 YEARS
(PHOTO BOTTOM LEFT) AND HIS 8 YEARS (PHOTO BOTTOM RIGHT).

**NORMALIZATION OF FUNCTIONS RESULTED IN A RESUMPTION OF EUMORPHIC
GROWTH AND SATISFACTORY DENTAL AND OCCLUSAL ALIGNMENT.**



FROGGY MOUTH