

**ISOLATED OPEN BITE:
REHABILITATION WITH
FROGGYMOUTH**
Dr. Patrick Fellus



31st of March



16th of May

THIS FOUR-YEAR-OLD GIRL PRESENTED WITH AN ISOLATED INCISOR OPEN BITE.

SHE WORE THE DEVICE BETWEEN MARCH 31 AND JUNE 1.



20th of June



11th of October

THE DYSMORPHOSIS CONTINUED TO CORRECT SPONTANEOUSLY, AND THE RESULT OBTAINED ON OCTOBER 11 SUGGESTS A PHYSIOLOGICAL OCCLUSION IN PERMANENT TEETH WITHOUT ORTHODONTIC ADJUSTMENT.



FROGGYMOUTH