

WE
CARE FOR
ENAMEL





DENTA
SONIC

Packaging refills, set and Kit's for your orders







Packaging & Refills - Aligner, Orthodontic, Lingual

DiaStrip 2 sides and 1 side.




Name:	DiaStrip	DiaStrip	DiaStrip	DiaStrip	DiaStrip
PN: 2 sides:	DSM-01-12	DSM-02-12	DSM-03-12	DSM-04-12	RDS-05-12
PN: 1 side left:	DSM-01L-12	DSM-02L-12	DSM-03L-12		
PN: 1 side right:	DSM-01R-12	DSM-02R-12	DSM-03R-12		
Thickness mm:	0.1mm	0.2mm	0.3mm	0.4mm	0.5mm

Diamond Resin DiaStrip 2 sides and 1 side.



Name:	Resin DiaStrip	Resin DiaStrip	Resin DiaStrip	Resin DiaStrip	Resin DiaStrip
PN: 2 sides:	RDS-02-12	RDS-03-12	RDS-04-12	RDS-05-12	RDS-06-12
Thickness mm:	0.2mm	0.25mm	0.3mm	0.4mm	0.5mm
PN: 1 side left:	RDS-008L-12	RDS-009L-12	RDS-010L-12	RDS-012L-12	
PN: 1 side right:	RDS-008R-12	RDS-009R-12	RDS-010R-12	RDS-012R-12	
Thickness mm:	0.08mm	0.09mm	0.10mm	0.12mm	

DiaSaw 2 sides and 1 side.



Name:	DiaSaw	DiaSaw	DiaSaw
PN: 2 sides:	DSS-01-12	DSS-02-12	DSS-03-12
Thickness mm:	0.01mm	0.02mm	0.05mm
PN: 1 side left:	DSS-005L-12		not coated
PN: 1 side right:	DSS-005R-12		
Thickness mm:	0.05mm		

DiaSaw no Diamond.

DiaBee 2 sides.



Name:	DiaBee	DiaBee	DiaBee	DiaBee	DiaPlaque
PN: 2 sides:	DSB-15-12	DSB-25-12	DSB-40-12	DSB-60-12	DP-1-12
Grain µm:	15µm	25µm	40µm	60µm	not coated

DiaPlaque no Diamond.

To order blister of only 1 or 3 strips, change the end of the code, from -12 to -1 or -3

Packaging & Refills - Aligner, Orthodontic, Lingual

DiaStrip 2 sides and 1 side.



Name:	DiaStrip	DiaStrip	DiaStrip	DiaStrip	DiaStrip
PN: 2 sides:	DSM-01	DSM-02	DSM-03	DSM-04	RDS-05
PN: 1 side left:	DSM-01L	DSM-02L	DSM-03L		
PN: 1 side right:	DSM-01R	DSM-02R	DSM-03R		
Thickness mm:	0.1mm	0.2mm	0.3mm	0.4mm	0.5mm

Diamond Resin DiaStrip 2 sides and 1 side.



Name:	Resin DiaStrip	Resin DiaStrip	Resin DiaStrip	Resin DiaStrip	Resin DiaStrip
PN: 2 sides:	RDS-02	RDS-025	RDS-03	RDS-04	RDS-05
Thickness mm:	0.2mm	0.25mm	0.3mm	0.4mm	0.5mm
PN: 1 side left:	RDS-008L	RDS-009L	RDS-010L	RDS-012L	
PN: 1 side right:	RDS-008R	RDS-009R	RDS-010R	RDS-012R	
Thickness mm:	0.08mm	0.09mm	0.10mm	0.12mm	

DiaSaw 2 sides and 1 side.



Name:	DiaSaw	DiaSaw	DiaSaw
PN: 2 sides:	DSS-01	DSS-02	DSS-1
Thickness mm:	0.01mm	0.02mm	0.05mm
PN: 1 side left:	DSS-0.05L		not coated
PN: 1 side right:	DSS-0.05R		
Thickness mm:	0.05mm		

DiaSaw no Diamond.



DiaBee 2 sides.



Name:	DiaBee	DiaBee	DiaBee	DiaBee	DiaPlaque
PN: 2 sides:	DSB-15	DSB-25	DSB-40	DSB-60	DP-1
Grain µm:	15µm	25µm	40µm	60µm	not coated

PoliStrip 2 sides and 1 side. Final polishing



Name: PoliStrip
PN: 2 sides: RPS-9
PN: 1 side left: RPS-9L
PN: 1 side right: RPS-9R
Grain μm : 9 μm
Content: 12X RPS-9

Name: PoliStrip
PN: 2 sides: RPS-15
PN: 1 side left: RPS-15L
PN: 1 side right: RPS-15R
Grain μm : 15 μm
Content: 12X RPS-15


PoliStrip 2 sides and 1 side. Final polishing



Name: PoliStrip
PN: 2 sides: RPS-S2
Grain μm : 9 μm , 15 μm
Content: 6X RPS-9, 15

Name: PoliStrip
PN: 2 sides: RPS1-S2
PN: 1 side left: RPS1-S2L
PN: 1 side right: RPS1-S2R
Grain μm : 9 μm , 15 μm
Content: 6X L or R RPS-9, 15


Polishing enamel surface result



Untreated enamel surface and treated surface with Resin PoliStrip 15 μm SEM enlargement 500x



Untreated enamel surface and treated with Resin PoliStrip 15 μm SEM enlargement 2000x




Treated enamel surface (40 μm) and center treated with Resin PoliStrip 15 μm SEM enlargement 2000x



Courtesy of the Institute of Material Science:
 Technische Universität Darmstadt, Germany

How insert and remove correctly the strip from the handpiece

Wrong strip handling



Never push inside the strip into the head like the two picture on the left, because you will not have enough force, and you risk to break the shank /blade of the strip.

Correct strip handling



To insert the strip proceed by push the shank into the head, by turning on the right and left side until the grip are ok.

To remove the strip proceed by pull out from the head by turning on the right and left side.

Never plie the strip because you can broke the shank. The shank resist until 1kg force by plie it.

Packaging & Set's - Aligner, Orthodontic, Lingual

DiaStrip 2 sides and 1 side.



Name:	DiaStrip	DiaStrip	DiaStrip	DiaStrip
PN: 1/2 sides:	DSM-S12	DSM-S5	DSM1-S12	DSM1-S5
Content:	2X DSM-0.1, 0.2, 0.3, 0.4 2X RDS-0.5, 2X RPS-15	1X DSM-0.1, 0.2, 0.3, 0.4, 1X RDS-0.5	2X DSM-0.1L, 0.1R, 0.2L, 0.2R, 0.3L, 0.3R	1X L or R DSM-0.1, 0.2, 0.3 2X RPS-15

Resin DiaStrip 2 sides.



Name:	Resin DiaStrip
PN: 1/2 sides:	RDS-S12
Content:	2X RDS-0.2, 0.25, 0.3, 0.4, 0.5 1X RPS-9, 15

DiaSaw 2 sides and 1 side.



Name:	DiaSaw	DiaSaw	DiaBee
PN: 1/2 sides:	DSS-S12	DSS-S5	DSB-S5
Content:	1X DSS-0.05L, 0.05R, 4X 0.1, 2X 0.2, 2X RPS-9, 15	1X DSS-0.05L, 0.05R, 0.1, 0.2, 1X RPS-15	1X DSB-15, 25, 40, 60 1X RPS-15

DiaBee 2 sides.



Aligner DiaStrip 2 sides.



Name:	Aligner DiaStrip	Aligner DiaStrip
PN: 1/2 sides:	DSA-S12A	DSA-S5A
Content:	2X DSS-0.1, 2X DSM-0.2, 0.3 1X DSM-0.4, 0.5 1X RDS-0.4, 0.5 2X RPS-15	1X DSS-0.1, 1X DSM-0.2, 0.3, 0.4 1X RDS-0.5

ErgoStrip Packaging & refills - Aligner, Orthodontic, Lingual

ErgoStrip 2 sides and 1 side



Name:	ErgoStrip	ErgoStrip	ErgoStrip	ErgoStrip
PN: 2 sides:	ESM-15-12	ESM-25-12	ESM-40-12	ESM-60-12
Content:	12X ESM-15	12X ESM-25	12X ESM-40	12X ESM-60
PN: 1 side:	ESM1-15-12	ESM1-25-12	ESM1-40-12	ESM1-60-12
Content:	12X ESM1-15	12X ESM1-25	12X ESM1-40	12X ESM1-60
Grain µm:	15µm	25µm	40µm	60µm

ErgoSaw 2 sides.



Name:	ErgoSaw	ErgoSaw	ErgoSaw	ErgoSaw
PN: 2 sides:	ESS-15-12	ESS-25-12	ESS-40-12	ESS-1
Content:	12X ESS-15	12X ESS-25	12X ESS-40	12X ESS-1
Grain µm:	15µm	25µm	40µm	not coated

ErgoBee 2 sides.



Name:	ErgoBee	ErgoBee	ErgoBee	ErgoBee
PN: 2 sides:	ESB-15-12	ESB-25-12	ESB-40-12	ESB-60-12
Content:	12X ESB-15	12X ESB-25	12X ESB-40	12X ESB-60
Grain µm:	15µm	25µm	40µm	60µm

until end of stock

ErgoPlaque no Diamond.



Name:	ErgoPlaque
PN: 1/2 sides:	EP-1
Content:	12X EP-1
Grain µm:	not coated



Ergostrip Refill & Set confection.

ErgoStrip Packaging & Set's - Aligner, Orthodontic, Lingual

ErgoStrip 2 sides and 1 side.



Name:	ErgoStrip	ErgoStrip	ErgoStrip
PN: 1/2 sides:	ESM-S12	ESM1-S12	ESM-S21
Content:	3X ESM-15, 25, 40, 60	3X ESM1-15, 25, 40, 60	1X ESM1-15, 25, 40, 60 2X ESM-15, 25, 40, 60

ErgoSaw 2 sides.



Name:	ErgoSaw
PN: 1/2 sides:	ESS-S12
Content:	4X ESS-15, 25, 40

ErgoBee 2 sides.



Name:	ErgoBee
PN: 1/2 sides:	ESB-S12
Content:	3X ESB-15, 25, 40, 60 <i>until end of stock</i>

Packaging & Refills - Restoration, Veneers, Crown & Bridges

Polishing ErgoStrip 2 sides, and 1 side.



Name:	Polishing ErgoStrip
PN: 1/2 sides:	REPS-15
Content:	12X REPS-15
Grain µm:	15µm

How use ErgoStrip



Finger Grip design for better tactile control & perfect Grip.

ErgoStrip Packaging & Refills - Aligner, Orthodontic, Lingual

ErgoStrip 2 sides and 1 side



Name:	ErgoStrip	ErgoStrip	ErgoStrip	ErgoStrip
PN: 2 sides:	ESM-15-4	ESM-25-4	ESM-40-4	ESM-60-4
Content:	4X ESM-15	4X ESM-25	4X ESM-40	4X ESM-60
PN: 1 side:	ESM1-15-4	ESM1-25-4	ESM1-40-4	ESM1-60-4
Content:	4X ESM1-15	4X ESM1-25	4X ESM1-40	4X ESM1-60
Grain µm:	15µm	25µm	40µm	60µm

ErgoSaw 2 sides.



Name:	ErgoSaw	ErgoSaw	ErgoSaw	ErgoSaw
PN: 2 sides:	ESS-15-4	ESS-25-4	ESS-40-4	ESS-1-4
Content:	4X ESS-15	4X ESS-25	4X ESS-40	4X ESS-1
Grain µm:	15µm	25µm	40µm	not coated

ErgoBee 2 sides.



Name:	ErgoBee	ErgoBee	ErgoBee	ErgoBee
PN: 2 sides:	ESB-15-4	ESB-25-4	ESB-40-4	ESB-60-4
Content:	4X ESB-15	12X ESB-25	4X ESB-40	4X ESB-60
Grain µm:	15µm	25µm	40µm	60µm
	<i>until end of stock</i>			

ErgoStrip Packaging & Set's - Aligner, Orthodontic, Lingual

ErgoStrip 2 sides and 1 side.



Name:	ErgoStrip	ErgoStrip
PN: 1/2 sides:	ESM-S4	ESM1-S4
Content:	1X ESM-15, 25, 40, 60	1X ESM1-15, 25, 40, 60

ErgoSaw 2 sides.



Name:	ErgoSaw
PN: 1/2 sides:	ESS-S4
Content:	1X ESS-15, 25, 40, 1

ErgoBee 2 sides.



Name:	ErgoBee
PN: 1/2 sides:	ESB-S4
Content:	1X ESB-15, 25, 40, 60
	<i>until end of stock</i>

Basic Set's

Basic Set Ortho

Art. Ref.:	DSG4-204
Stroke:	4:1
Pieces:	14
Sterilizable:	Yes
Application:	■ ■
with a 1 additional reciprocating head	
Warranty option:	1 or 2 years



1 additional reciprocating head

Silent handpiece, no vibration. Internal water cooling and interchangeable head

Basic Set Aligner

Art. Ref.:	DSG4A-204
Stroke:	4:1
Pieces:	11
Sterilizable:	Yes
Application:	■ ■
Warranty option:	1 or 2 years



Basic Set Orthofamily

Art. Ref.:	DSG4-204-OF
Stroke:	4:1
Pieces:	13
Sterilizable:	Yes
Application:	■ ■
Warranty option:	1 or 2 years



Tested by Tor Vergata University, Rome.

Recognised Invisalign center for Invisalign certification courses.

Basic Set Ortho NWC

Art. Ref.:	DSG4-200
Stroke:	4:1
Pieces:	12
Sterilizable:	Yes
Application:	■ ■
Warranty option:	1 year



No Water-cooled Handpiece

Professional Kit 1

Art. Ref.:	DSPK1-2Y
Stroke:	4:1
Pieces:	35
Sterilizable:	Yes

with 1 additional rotating push button head
Warranty option: 2 years



Professional Kit 2



Art. Ref.:	DSPK2-2Y	DSPK2-1Y	DSPK2L-2Y
Stroke:	4:1	4:1	4:1
Pieces:	35	35	35
Sterilizable:	Yes	Yes	Yes

with 1 additional reciprocating head,
and voucher

Warranty option: 2 years 1 year 2 years



Professional Kit 3

Art. Ref.:	DSPK3-2Y
Stroke:	4:1
Pieces:	36
Sterilizable:	Yes

with 1 additional reciprocating and 1 rotating head
Warranty option: 2 years





DENTA
SONIC

Take
the guesswork
out of

IPC interproximal
correction



DENTASONIC

For more information about all this products, please contact us online:
5% off for all online order with code "DS5" on www.denta-sonic.com

DS DENTASONIC SA
Via Piodella 12
CH-6933 Muzzano
Switzerland

Tel. +41 91 600 19 60
info@dentasonic.ch
www.dentasonic.ch

Official Distributor:

European Representative: QuaSer S.r.l., Via Pian due Torri, 19, 00146 Roma

Some of the products and designations mentioned in the text are patented or copyrighted. The absence of a special reference or the ® sign should not be interpreted as the absence of copyright protection. This publication is copyrighted. All rights, including those of translation, reprinting, and reproduction (also in the form of extracts) are reserved. No part of this publication may be reproduced or processed using electronic systems in any form or by any means (photocopying, microfilm, or other methods) without the editor's prior consent in writing. Colours and products subject to alterations. Printing errors excepted.